

# Imprimis Wipe

Disinfecting/ Deodorizing/  
Cleaning

### Highlights

A non-woven disposable cloth containing a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces where control of cross-contamination is required.

This product is effective against the following pathogens in 2 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed:

Mycobacterium tuberculosis, Human Hepatitis B virus, Human Hepatitis C virus, Herpes simplex virus 2 (genital Herpes virus), HIV-1 (AIDS virus)\* (At an exposure time of 1 min.), \*Influenza A (A2 Japan), Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella enterica, Methicillin Resistant Staphylococcus aureus (MRSA), Escherichia coli O157:H7, Vancomycin Resistant Enterococcus faecalis (VRE).

**PREPARATION OF DISPENSER:** Remove cover from canister and discard covering seal. From center of wipes roll, pull up a wipe corner, twist into a point and thread through the X hole in the underside of the lid cover. Pull wipe through about one (1) inch. Replace the lid cover on container, assuring a secure fit. Pull out first wipe and snap at a 90° angle. Remaining wipes feed automatically for next use. Always reseal hinged cap after each use to insure moisture retention.

**DISINFECTION:** Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 2 minutes and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, nonporous, inanimate surfaces such as stainless steel, Formica, glass tables, chairs, baskets, counter tops, exam tables, toilets, basins, tubs, bedside, chart racks, equipment, cabinets, door knobs, sinks, mirrors, light fixtures, carts, waste cans, laundry hampers, wheel chairs, walkers, urinals, plastic and other similar hard surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough rewiping with a fresh wipe is required for disinfection. Use on food contact surfaces, must be followed by a potable water rinse.

**DEODORIZING:** Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, fecal matter, urine, burns, smoke, mildew and mold, decay, stale air, tobacco, sickrooms and musty closets.

**CLEANING:** To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. Use in bathrooms, sickrooms. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

### Applications

Hospitals, intensive care units, surgery, recovery, anesthesia, X-ray, Cat. lab., orthopedics, newborn nursery, respiratory therapy, emergency medical settings, laboratories, clinics, nursing facilities, medical offices.

### Active Ingredients

N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)	
Dimethyl Benzyl Ammonium Chloride .....	0.12%
N-Alkyl (68% C12, 32% C14)	
Dimethyl Ethyl Benzyl Ammonium Chloride .....	0.12%
Isopropyl Alcohol .....	41.58%
Inert Ingredients† .....	58.18%
TOTAL.....	100.00%

†Does not include weight of dry wipe.