

Certus Medical brings years of product development experience to design a line of products that are intended for process improvement. These solutions will provide for an exceptional reprocessing experience when clean instruments are utmost importance.

Not every Sterile Processing Department is the same. Certus Medical can help to design the right group of products that meet your facility's needs.

Virtus Enzymatic Transport Foam and Gel

Improved outcomes begin with instruments that are easy to clean once arriving in to the sterile processing department. Virtus Enzyme Transport Spray and Foam are designed to make that happen. These dual enzyme transport spray and foam products are uniquely formulated to keep the instruments wet from the moment they are placed in the tray to the time they arrive in the reprocessing unit. The multiple enzymes, including Lipase, are essential to break down the fat and other bioburden on the instruments, making it easier for Central Processing personnel to manage instrumentation. This is a crucial first step in reprocessing instrumentation.

Available Package Sizes

12/24 ounce case	697424PM
12/24 ounce case	697424PMA

Virtus Enzymatic Cleaner

Contains multiple enzymes, including Lipase. This product works exceptionally well in a broad range of water conditions and as an effective step in the breakdown of bioburden and debris. Highly effective cleaner perfect for automated washers or manual use.

Available Package Sizes

4/1 gallon case	662241PM
5 gallon pail	662205PM
15 gallon drum	662215PM

Virtus Plus Fast Acting Multi-enzyme

High levels of three enzymes are formulated for the most challenging cases. Effective at removing fat, blood, and other debris from all types of instruments. Works in all ranges of water conditions and is designed to be an effective step in the reprocessing of instrumentation. Safe for automated washers and for manual use.

Available Package Sizes

4/1 gallon case	660341PM
5 gallon pail	660305PM
15 gallon drum	660315PM

Integra Neutral pH Detergent

A broad based neutral pH detergent that is useful in a full range of water conditions. Use along with the Certus brand enzymes to provide an effective reprocessing program. Neutral pH is useful for automated and manual reprocessing.

Available Package Sizes

4/1 gallon case	697341PM
5 gallon pail	697305PM
15 gallon drum	697315PM

Integra Instrument Lubricant

Integra Instrument Lubricant is a valuable part in the reprocessing of instrumentation. Formulated to eliminate stiffness at the box locks of scissors and other hinged instruments. Can be used as a manual spray on the clean side of the sterile processing area or as a concentrate for use in automated reprocessing machines. Available in ready to use or as a concentrate.

Available Package Sizes

4/1 gallon case	697241PM
5 gallon pail	697205PM
15 gallon drum	697215PM
12/24 ounce case	661624PM

Eluo Spot Free Rinse

In many geographical areas the water is not always conducive for easy cleaning of instrumentation. Eluo Spot Free Rinse Dry is an important final step. Eluo is non-sudsing. Spotless drying without wiping! For use on glass, plastics, and all types of instrumentation. Eluo is a concentrated liquid rinse additive for use with injection pumps or in manual wash/rinse operations.

Available Package Sizes

5 gallon pail	618805PM
---------------	----------

ACI 100

Formulated with a blend of phosphoric and Hydroxyacetic acids and non-foaming detergent for use in automated reprocessors. Removes hard water scale, rust and other mineral deposits. Brightens stainless steel and aluminum equipment. Neutralizes alkali residues when used for acidified rinse. Useful to clean the inside of automated washers.

Available Package Sizes

5 gallon pail	651005PM
---------------	----------

Alka 150

Non-foaming, safe for all types of instrumentation including soft metals. This concentrated non-foaming cleaner detergent is designed for use in automated machines for removal of debris and other soils from bed pans, utensils, carts, and racks. Useful as a cart washer.

Available Package Sizes

5 gallon pail	653805PM
---------------	----------

Stainless Steel Maintainer

A water based cleaner used for polishing and cleaning stainless steel, chrome aluminum, plastics and other hard surfaces. Does not have the oiliness and flammability of solvent based stainless steel cleaners. Solvents have been removed and replaced with water to lower the VOC's and make it more environmentally friendly.

Available Package Sizes

12/20 ounce case	458720PM
------------------	----------

