



6622

Virtus Enzymatic Cleaner

Ingredients

Intentionally Added Ingredients

<u>CAS number</u>	<u>Chemical name</u>	<u>Role</u>	<u>Designated Lists</u>
7732-18-5	Water	solvent	
67-63-0	Isopropanol	solvent	CA Non-Cancer Hazard
126-92-1	Sodium Ethylhexyl Sulfate	surfactant	
not available	Sodium Alkanoate	surfactant	
102-71-6	Triethanolamine	pH adjuster	
9014-01-1	Protease Enzyme	Biocatalyst	EU Respiratory Sensitizer
64-18-6	Formic Acid	stabilizer	
9001-62-1	Lipase Enzyme	Biocatalyst	
64665-57-2	Sodium Tolyltriazole	corrosion inhibitor	
10043-52-4	Calcium Chloride	stabilizer	
3844-45-9	FD&C Blue 1	colorant	
633-96-5	D&C Orange 4	colorant	
4403-90-1	D&C Green 5	colorant	
3844-45-9	Acid Blue 9	colorant	

Nonfunctional Constituents

<u>CAS number</u>	<u>Chemical name</u>	<u>Role</u>	<u>Designated Lists</u>
	none		