

Integra

Instrument Lubricant
Concentrate

INTEGRA LUBRICANT Integra Instrument Lubricant is a valuable part in the reprocessing of instrumentation. Formulated to eliminate stiffness at the box locks of scissors and other hinged instruments. Can be used as a manual spray on the clean side of the sterile processing are or as a concentrate for use in automated reprocessing machines.

Highlights

- Increase instrument life
- Keeps scissors sharp
- Safe on all types of metals
- Reduces rusting, corrosion and discoloration.

Typical Properties

- Specific Gravity: 0.95
- pH: 7.5
- Appearance: Milky white emulsion
- Fragrance: odorless

Applications

For cleaning and lubricating surgical instruments and medical devices

Safety Cautions

CAUTION: Eye and skin irritant. Keep out of reach of children. Refer to M.S.D.S. for additional safety information.