

Virtus Enzymatic Transport Gel

Pre-Cleaner Spray Application

Virtus Enzymatic Transport Gel is formulated to keep the instruments wet from the moment they leave each department to the time they arrive in the reprocessing unit. The multiple enzymes, including Lipase, are essential to break down the fat and other bioburden on the instruments, making it easier for Central Processing personnel to manage them. This is a crucial first step in reprocessing instrumentation.

Highlights

- Multiple Enzymes formula
- Retains moisture
- Contains corrosion inhibitors
- Minimizes exposure to soil

Typical Properties

- Specific Gravity: 1.0
- pH: 8.0
- Appearance: Clear liquid
- Fragrance: Odorless

Applications

- Hospitals, nursing homes (patient turnover, physical therapy rooms),
- Hospice care centers
- Doctors offices
- Surgery centers
- Acute care centers

Safety Cautions

Keep out of reach of children. Refer to M.S.D.S. for additional safety information.